



Floor Safety & Protection




Xodus Medical proudly supports the **AORN Tool Kit** for the Prevention of Perioperative Pressure Injury. Scan to view up-to-date patient safety recommendations or visit xodusmedical.com/AORNtoolkit



Camel Mats | Absorbent Floor Protection

Solutions for Maintaining Clean, Dry Hospital Floors




- Keep hospital floors clean and dry
- Single-use design improves infection control and reduces cross-contamination risk
- Save money - reduce turnover time
- Available in different sizes and absorbency levels
- Non-slip backing for secure placement



Non-Drip Disposal

The Mojave Max mat features a polymer material that enhances absorption and prevents fluid



Advanced Fluid Retention Technology

Fluid-absorbing polymer traps liquid for easy and convenient clean up.

Mojave Standard Absorbent Mat



Prod #	Size
80510	24" x 32"

Mojave Ultra Absorbent Mat



Prod #	Size
80512	24" x 32"
80522	32" x 40"
80526	32" x 56"
80524	48" x 80"
80525	48" x 80"

-w/Cutout

Mojave Max Absorbent Mat



Prod #	Size
80532	24" x 32"
80531	32" x 40"
80530	32" x 56"

Camel® Cord Cover

Two peel strips on the back reveal residue-free adhesive that secures the cover to the floor. This feature safely secures cords and decreases the risk of trips and falls. Available in two sizes.




#80548



Soft Step™ | Anti-Fatigue Floor Mats

Increasing Safety and Comfort for Hospital Staff

- Tapered channels prevent problematic puddles in the OR
- Decrease muscle fatigue
- Increase safety - keep clinical staff comfortable
- Available in extra-long and mini sizes



Product #	Size
80543	70" x 20"
80544	35" x 20"
80545	17" x 13"

#80545



#80541



#80537

"Prolonged standing during surgical procedures poses a high risk of causing musculoskeletal disorders, including back, leg, and foot pain, which can be chronic or acute in nature."

— AORN Journal. June, 2011

Fitted Camel Mat

Extend the usable life of Soft Step mats by incorporating single-use, absorbent covers.

- Available in two sizes:
#80536 - 39" x 24"
#80537 - 74" x 24"




The Oasis™ | Anti-Fatigue Floor Mats

Increasing Safety and Comfort for Hospital Staff



The Oasis™ Mat deters the onset of muscle fatigue and helps keep hospital floors clean, dry, and safe.

At the core of the Oasis, between two layers of shock absorbing, closed-cell foam, is the same ultra-absorbent material featured in the Camel™ Mat. An impermeable barrier prevents fluid from reaching the hospital floor, while a non-slip backing keeps the mat securely in place.



- Maximize performance by decreasing muscle fatigue
- Improve infection control by reducing exposure to contaminants
- Increase safety by minimizing possibility of injury from slips and falls
- Save money by reducing turnover time

Stacking up the Protection



XODUS MEDICAL
Making Surgery Safer™

702 Prominence Drive
New Kensington, PA 15068,
United States

xodusmedical.com
info@xodusmedical.com
Tel: 1+724-337-5500

FDA REGISTERED
ISO 13485 CERTIFIED

LITBRAM006 R4 1/11/2024



MADE IN THE
USA



Follow Us on
Social Media