




Floor Safety & Protection



Camel Mats™ | Absorbent Floor Protection

Solutions for Maintaining Clean, Dry Hospital Floors




Help Eliminate Slips, Trips, and Falls

- Keep hospital floors clean and dry
- Single-use design improves infection control and reduces cross-contamination
- Non-slip backing for secure placement

Advanced Fluid Retention Technology

Fluid-absorbing polymer traps liquid for easy and convenient cleanup and enhances absorption, preventing fluid drips



Camel Cord Cover™

Two peel strips reveal residue-free adhesive that secures the cover to the floor. This feature safely secures cords and decreases the risk of trips and falls.

Available in two sizes



#80548

Camel Mat – Standard Absorbent



Prod

#80510
#80520
#80527

Size

24" x 32"
32" x 40"
32" x 56"

Camel Mat – *Ultra* Absorbent



Prod

#80512
#80522
#80526
#80524
#80525

Size

24" x 32"
32" x 40"
32" x 56"
48" x 80"
*48" x 80"
*w/ cutout

Camel Mat – *Max* Absorbent




Prod

#80532
#80531
#80530

Size

24" x 32"
32" x 40"
32" x 56"



#80525
w/ cutout

Soft Step™ | Anti-Fatigue Floor Mats


Increasing Safety and Comfort for Hospital Staff

- Tapered channels prevent problematic puddles in the OR
- Decrease muscle fatigue
- Increase safety - keep clinical staff comfortable
- Available in extra-long and mini sizes




#80545

Product #	Size
80543	70" x 20"
80544	35" x 20"
80545	17" x 13"



#80537



#80541

Fitted Camel Mat

Extend the usable life of Soft Step mats by incorporating single-use, absorbent covers

- Available in two sizes:
#80536 - 39" x 24"
#80537 - 74" x 24"

// "Prolonged standing during surgical procedures poses a high risk of causing musculoskeletal disorders, including back, leg, and foot pain, which can be chronic or acute in nature."


— AORN Journal. June, 2011

The Oasis™ | Anti-Fatigue Floor Mats

Increasing Safety and Comfort for Hospital Staff



The Oasis™ Mat deters the onset of muscle fatigue and helps keep hospital floors clean, dry, and safe.

At the core of the Oasis, between two layers of shock absorbing, closed-cell foam, is the same ultra-absorbent material featured in the Camel™ Mat. An impermeable barrier prevents fluid from reaching the hospital floor, while a non-slip backing keeps the mat securely in place.



- Maximize performance by decreasing muscle fatigue
- Improve infection control by reducing exposure to contaminants
- Increase safety by minimizing possibility of injury from slips and falls
- Save money by reducing turnover time

Stacking up the Protection



XODUS MEDICAL
Making Surgery Safer™

702 Prominence Drive
New Kensington, PA 15068,
United States

xodusmedical.com
info@xodusmedical.com
Tel: 1+724-337-5500

FDA REGISTERED
ISO 13485 CERTIFIED

LITBRAM006 R7 2/18/2025



MADE IN THE
USA



Follow Us on
Social Media