



AORN

PATIENT SAFETY

EFFICIENCY  
QUALITY  
SAFETY

SAFER SURGERY

ADVANCED

SUPERIOR PROTECTION

QUALITY

PROPER POSITIONING

MINIMIZE RISK

PREVENTION  
INNOVATION

# Sterile Processing Solutions




Xodus Medical proudly supports the **AORN Tool Kit** for the Prevention of Perioperative Pressure Injury. Scan to view up-to-date patient safety recommendations or visit [xodusmedical.com/AORNtoolkit](http://xodusmedical.com/AORNtoolkit)



# GreenGuard Tray Protectors

Reusable Solutions for Surgical Tray Protection



## Simplify Tray Protection

- Ideal for heavy trays, both steam and Sterrad compatible, works equally well on square or round-corner trays

## Reduce Reprocessing Time

- Greater airflow decreases drying time and eliminates moisture retention

## Increase Patient Safety

- Imperfections in sterile wrapping provide opportunities for contamination, thereby increasing the risk of SSI<sup>1</sup>. Because sterile wrap perforations are often overlooked in the Operating Room, prevention of tears or perforations is critical<sup>2</sup>

## Protect the Environment

- Reusable, recyclable, and reduces the need for extra sterile wraps

**“Minimum re-sterilization “wait cost” for the OR [is] \$1,342 (per incident)... Instrument availability and reliability are essential to a well-run, cost-effective OR environment.”<sup>3</sup>**

1. Looking for Holes in Sterile Wrapping: How Accurate Are We? <https://doi.org/10.1007/s11999-0000000000000185>  
2. Implant contamination as a cause of surgical site infection in spinal surgery: are single-use implants a reasonable solution? - A systematic review. <https://doi.org/10.1186/s12891-020-03653-z>  
3. The Clinical and Economic Case for Sterile, Disposable Instruments and Implants. <https://www.infectioncontroltoday.com/view/clinical-and-economic-case-sterile-disposable-instruments-and-implants>



# Single-Use Corner Protectors


Single-Use Solutions for Surgical Tray Protection



#40461

## Prevent Tears & Moisture In Sterile Wraps

- Create a soft barrier to protect sterilization wrap from sharp corners; prevent holes and/or tears
- Eliminate unnecessary costs created by contaminated surgical trays (*Average savings: \$40 per contaminated tray*)
- Avoid double-wrapping trays
- Cost-effective solution




#40463

**40461** White foam 1.25" thick  
100 sheets/case (4 squares | 6" x 6" sheet)

**40462** Reticulated charcoal foam 1.25" thick  
100 sheets/case (4 squares | 6" x 6" sheet)

**40463** Reticulated charcoal foam 1.625" thick  
80 sheets/case (4 squares | 6" x 6" sheet)

## Depth Comparison





# Camel Mats | Absorbent Floor Protection

Solutions for Maintaining Clean, Dry Hospital Floors



**Absorbent Mat**

**Fluid-Proof Poly Barrier**  
(non-slip versions)

**Non-Slip Backing**

- Keep hospital floors clean and dry
- Single-use design improves infection control and reduces cross-contamination risk
- Save money - reduce turnover time
- Available in different sizes and absorbency levels
- Non-slip backing for secure placement

## Mojave Standard Absorbent Mat



**Prod #**  
80510

**Size**  
24" x 32"

## Mojave Ultra Absorbent Mat



**Prod #**  
80512  
80522  
80526  
80524  
80525

**Size**  
24" x 32"  
32" x 40"  
32" x 56"  
48" x 80"  
48" x 80"  
-w/Cutout

## Mojave Max Absorbent Mat



**Prod #**  
80532  
80531  
80530

**Size**  
24" x 32"  
32" x 40"  
32" x 56"




#80548



### Non-Drip Disposal

The Mojave Max mat features a polymer material that enhances absorption and prevents fluid from dripping.



### Advanced Fluid Retention Technology

Fluid-absorbing polymer traps liquid for easy and convenient clean up.



#80548

### Camel® Cord Cover

Two peel strips on the back reveal residue-free adhesive that secures the cover to the floor. This feature safely secures cords and decreases the risk of trips and falls. Available in two sizes.




#80512



# Soft Step™ | Anti-Fatigue Floor Mats

Increasing Safety and Comfort for Hospital Staff

- Tapered channels prevent problematic puddles in the OR
- Decrease muscle fatigue
- Increase safety - keep clinical staff comfortable
- Available in extra-long and mini sizes



Product #	Size
80543	70" x 20"
80544	35" x 20"
80545	17" x 13"

#80545



#80541



#80537

“Prolonged standing during surgical procedures poses a high risk of causing musculoskeletal disorders, including back, leg, and foot pain, which can be chronic or acute in nature.”

— AORN Journal. June, 2011

## Fitted Camel Mat

Extend the usable life of Soft Step mats by incorporating single-use, absorbent covers.

- Available in two sizes:  
#80536 - 39" x 24"  
#80537 - 74" x 24"




# The Oasis™ | Anti-Fatigue Floor Mats

## Increasing Safety and Comfort for Hospital Staff



The Oasis™ Mat deters the onset of muscle fatigue and helps keep hospital floors clean, dry, and safe.

At the core of the Oasis, between two layers of shock absorbing, closed-cell foam, is the same ultra-absorbent material featured in the Camel™ Mat. An impermeable barrier prevents fluid from reaching the hospital floor, while a non-slip backing keeps the mat securely in place.



- Maximize performance by decreasing muscle fatigue
- Improve infection control by reducing exposure to contaminants
- Increase safety by minimizing possibility of injury from slips and falls
- Save money by reducing turnover time

### Stacking up the Protection



**XODUS MEDICAL**  
Making Surgery Safer™

702 Prominence Drive  
New Kensington, PA 15068,  
United States

[xodusmedical.com](http://xodusmedical.com)  
[info@xodusmedical.com](mailto:info@xodusmedical.com)  
Tel: 1+724-337-5500

FDA REGISTERED  
ISO 13485 CERTIFIED



MADE IN THE  
**USA**



Follow Us on  
**Social Media**