





# Electrosurgical

Pencils, Electrodes,  
and Accessories





A Complete Line  
of Superior Quality  
Electrosurgical  
Products




# Laparoscopic Tips


## Product Numbers


	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20080		#20145
	Foot Control	#20127		#20111


	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20081		#20138
	Foot Control			

	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20082	#20157	#20146
	Foot Control	#20115		#20147

	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20083		#20148
	Foot Control	#20118		#20149

	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20106	#20158	#20150
	Foot Control			

	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20135		#20151
	Foot Control			














	Size	33 cm.	37 cm.	44 cm.
	Hand Control	#20144		#20152
	Foot Control	#20124		#20153

## Patient Safety First 100% Tested



- Industry leading, heat-resistant Kynar insulation
  - » Single-use eliminates reprocessing
  - » Eliminates cross-contamination
  - » Lowers total costs
  - » Prevents insulation damage, associated burn risk
- Compromised insulation may lead to stray currents and thermal necrosis of non-targeted tissue— which carries the potential for peritonitis, other serious injury or even death
- Every electrode individually tested (high-potential stress test)
- PTFE coatings deliver consistent coagulation and resist flaking at high temperatures



# PTFE Coated Electrodes

	6.5" PTFE Coated Needle Tip Blade Electrode	Product No: #20023	50 /Case
	6.5" Standard Blade Electrode	Product No: #20062	50 /Case
	6.5" PTFE Coated Standard Blade Electrode	Product No: #20014	50 /Case
	6.5" PTFE Coated Standard Blade Electrode	Product No: #20015	50 /Case
	5" PTFE Coated Ball Tip Blade Electrode	Product No: #20018	50 /Case
	4" PTFE Coated Standard Blade Electrode	Product No: #20012	50 /Case
	4" PTFE Coated Standard Blade Electrode	Product No: #20013	50 /Case
	2.75" PTFE Coated Standard Blade Electrode	Product No: #20057	50 /Case
	2.75" PTFE Coated Standard Blade Electrode	Product No: #20058	50 /Case
	2.75" PTFE Coated Standard Blade Electrode	Product No: #20055	50 /Case
	2.75" PTFE Coated Standard Blade Electrode	Product No: #20056	50 /Case
	2.5" PTFE Coated Standard Blade Electrode	Product No: #20010	50 /Case
	2.5" PTFE Coated Standard Blade Electrode	Product No: #20011	50 /Case
	6.25" PTFE Coated Bayonet Blade Electrode	Product No: #20017	50 /Case

# Tungsten Micro-Needles

	1.6" Tungsten Micro-Needle Blade Electrode	Product No: #20072	12 /Case
	2.75" Tungsten Micro-Needle Blade Electrode	Product No: #20070	12 /Case
	2.75" Tungsten Needle Tip Blade Electrode	Product No: #20102	12 /Case

# Pencils with PTFE Coated Electrodes


High quality components and construction allow for precise switch control, safeguarding against accidental activation. Performance is exceptional at lower power settings, helping to minimize thermal necrosis risk. Patient safety is maximized through strict quality control and the use of premium materials. Optional preloaded, PTFE coated electrodes minimize eschar build-up.



All pencils are available with PTFE or Stainless Tips  
View a full list of pencil and electrode combinations on [xodusmedical.com](http://xodusmedical.com)

Holsters Included With All Electrocautery Pencils


## Grounding Pads



- Round, latex-free, split design with and without a 10 ft (3m) cord
- Overlapping gel-line border limits aluminum edge exposure, reducing burn risk
- Compatible with generators featuring contact quality monitors

Product No.	Description	Case QTY
#20240	Split Grounding Pad with Cord	50
#20230	Split Grounding Pad without Cord	50

## Cautery Tip Cleaners



- Abrasive surface effectively removes tissue and debris with one pass
- Cleans uncoated (stainless steel) cauter tips
- Adhesive backing for universal placement
- X-ray detectable

Product No.	Case QTY
#30500	50

## EZ Covers



- Prevents cross-contamination
- Easy open/close draw string
- Easy cleanup of blood & other infectious materials

Product No.	Case QTY
#50300	50



**XODUS MEDICAL**  
Making Surgery Safer™

702 Prominence Drive  
New Kensington, PA 15068,  
United States

[xodusmedical.com](http://xodusmedical.com)  
[info@xodusmedical.com](mailto:info@xodusmedical.com)  
Tel: 1+724-337-5500

FDA REGISTERED  
ISO 13485 CERTIFIED



MADE IN THE  
**USA**



Follow Us on  
**Social Media**