




The Pink Lat Pack™

A New Standard for Safe, Effective, and Efficient Lateral Positioning


Multi-Angle Wedge

35°, 55°, and 90° inclination options ensure ideal positioning for every patient and procedure



Adjustable Lateral Head Positioner


Skin-friendly materials, adjustable height, and an ergonomic design make this the ideal Lateral Head Positioner



Breathable, Skin-friendly, Moisture-wicking Material

Safety Maximized

Increase patient safety with DermaProx™ and PinkProtect™: breathable, moisture-wicking materials that gently immerse and envelop areas prone to pressure injury.



Standard of Care Emphasized

No customization required. Pink Lat Packs™, regardless of weight or BMI, provide a safe, consistent standard of care for every patient, in every procedure, every time.



Efficiency Optimized

In clinical trials, an average of 15-20 minutes were saved during setup. When necessary, The Pink Lat Pack™ makes repositioning faster, safer, and easier.


“HAPI occurrence affects 2.5 million patients per year, adding more than \$43,180 per patient on average in-hospital costs.”

(October 2014. Agency for Healthcare Research & Quality)

“For years prior to the introduction of Pink Lat Packs, we at Yale New Haven Hospital relied on reusable devices for patient positioning in the lateral-decubitus position. Interrupting procedures to reposition patients could negatively impact both patient safety and surgical efficiency. Pink Lat Packs are now an important component of our efforts to safely secure our patients, protect their skin and, ultimately, improve surgical outcomes.”

— Dr. Domenico Savatta, MD
 Chief of Minimally Invasive & Robotic Adult Urologic Surgery
 New Jersey


Pink Lat Pack™ Available Configurations



*Not all components are pictured

#LPK10022 – Savatta Pink Lat Pack™


- A. Adjustable Lateral Head Positioner (1)
- B. Multi-Angle Body Wedge (1)
- C. Skin-Friendly Body Straps (4)
- D. Lateral Leg Positioner (Savatta Style) (1)
- E. Pink Pad EXT™ (1)
- F. Pink Arc (Axillary Roll) (1)
- G. One-Step™ Arm Protectors with DermaProx™ (1pr)
- H. Lift Sheet (1)



*Not all components are pictured

#LPK10005 – Standard Pink Lat Pack™


- A. Adjustable Lateral Head Positioner (1)
- B. Multi-Angle Body Wedge (1)
- C. Tapered Intra-Leg Positioner (1)
- D. Pink Pad EXT™ (1)
- E. Skin-Friendly Body Straps (4)
- F. Armboard Pads with DermaProx™ (1pr)
- G. Skin-Friendly Armboard Straps (1pr)
- H. Lift Sheet (1)



*Not all components are pictured

#LPK10023 – Jackknife Pink Lat Pack™

- A. Adjustable Lateral Head Positioner (1)
- B. Multi-Angle Body Wedge (2)
- C. Skin-Friendly Body Straps (3)
- D. Tapered Intra-Leg Positioner (1)
- E. Pink Pad EXT™ (1)
- F. Pink Arc (Axillary Roll) (1)
- G. Armboard Pads with DermaProx™ (1pr)
- H. Skin-Friendly Armboard Straps (1pr)
- I. Lift Sheet (1)




*Not all components are pictured

#LPK10015 – Thoracic Pink Lat Pack™

- A. Adjustable Lateral Head Positioner (1)
- B. Low Profile Lateral Wedge (2)
- C. Intra-Leg Positioner (1)
- D. Armboard Pads with PinkProtect™ (1pr)
- E. Skin-Friendly Arm Board Straps (2)
- F. Pink Pad EXT™ (1)
- G. Lift Sheet (1)
- H. Pink Arc (Axillary Roll) (1)
- I. Skin-Friendly Body Straps (3)


Individual Positioning Solutions



40941 - Lateral Leg Positioner (Savatta Style)



40471 - Adjustable Lateral Head Positioner



40600 - Pink Pad™ EXT (w/ Lift Sheet + Body Strap)




40943 - Small Body Wedge



40948 - Multi-Angle Body Wedge (Small)




40947 - Low-Profile Lateral Wedge




40945 - 3" Intra-Leg Positioner



40547 - Pink Armboard Pads w/ PinkProtect™ + Arm Straps



40946 - Tapered Intra-Leg Positioner



Various Options Available
PinkProtect™ Dual Layer Pads



Various Options Available
Pink Arcs



XODUS MEDICAL
Making Surgery Safer™

702 Prominence Drive
New Kensington, PA 15068,
United States

xodusmedical.com
info@xodusmedical.com
Tel: 1+724-337-5500

FDA REGISTERED
ISO 13485 CERTIFIED

LITBRLP001 R12 10/13/2023



MADE IN THE
USA



Follow Us on
Social Media