

# The Pink Reverse-T Kit™

Reverse Trendelenburg Positioning System



- **Safely Prevent Sliding**

An innovative seat contour made of high-density, supportive material ensures anatomically correct positioning, while skin-friendly PinkProtect™ safely prevents unwanted movement in reverse Trendelenburg

- **Reduce Risk of HAPIs**


Beyond preventing unwanted movement in reverse Trendelenburg, the Pink Reverse-T Kit features a unique material that offers patients unparalleled skin protection.

Unique softness, breathability, immersion, and envelopment characteristics protect against Hospital-Acquired Pressure Injuries (HAPIs).

- **Complete Coverage for Every Patient**


The Pink Reverse-T Kit™ requires less than one minute to set up. A single-use design minimizes cross-contamination risk and improves OR turnover times.

## Components




**#40519**  
PinkProtect™ Head  
Positioner

Multiple Pink Pad  
Configurations Available



**#40433**  
One-Step Arm Protectors  
(Optional)




**#40548**  
Reverse-T Footboard Pad  
(Optional)

#### U.S. Patents

8,464,720; 8,511,314; 9,161,876; 9,750,656; 9,782,287;  
9,931,262; 9,949,883; 10,045,902; 10,285,890; 10,324,050;  
2,859,398; 11,484,456


**Canada Patents:** 2,859,398  
Additional Patents Pending






“With patients weighing 450lbs or more, and BMIs in excess of 65, safe positioning can be a real challenge. Xodus Medical’s **Pink** Reverse-T Kit not only prevented our patients from sliding, but did so safely.”

— Alicia C, RN, BSN, Surgical Services Manager, Hartford, KY



### Popular Configuration (#40640) Includes:

- 1 Pink Reverse-T Pad  
w/ Posterior Support Contour
- Large Armboards (Pair) – with DermaProx
- 1 Lift Sheet
- 2 Sets of Skin-Friendly Body Straps
- Single Use Adult Head Positioner
- Reverse-T Footboard Pad (20"x12"x2")



### Popular Configuration (#40643) Includes:

- 1 Pink Pad XL
- Large Armboards (Pair) – with DermaProx
- 1 Lift Sheet
- 2 Sets of Skin-Friendly Body Straps
- Single Use Adult Head Positioner
- Reverse-T Footboard Pad (20"x12"x2")



**XODUS MEDICAL**  
Making Surgery Safer™

702 Prominence Drive  
New Kensington, PA 15068,  
United States

xodusmedical.com  
info@xodusmedical.com  
Tel: 1+724-337-5500

FDA REGISTERED  
ISO 13485 CERTIFIED



MADE IN THE  
**USA**



Follow Us on  
**Social Media**