

Sterile Processing Products



The Soft Step XL

- Tapered channels prevent problematic fluid pooling
- Decreases muscle fatigue
- Increases safety - keeps staff comfortable

PRODUCT	SIZE
#80543	70" X 20"
#80544	35" X 20"

"(Prolonged standing) poses a high risk of causing musculo-skeletal disorders, including back, leg, and foot pain, which can be chronic or acute in nature." (AORN Journal, June 2011)


FITTED CAMEL MAT

The Soft Step™ XL is double the length of the regular Soft Step™ to help prevent trips and falls from "doubling up" floor mats.

Available in two sizes:
 #80536 - 39" x 24"
 #80537 - 74" x 24"



A line of absorbent mats that feature multiple layers designed to protect both hospital floors and staff. Available in three different levels of absorbency.



The Gobi



The Mojave




The Sahara

Camel Mats come in a variety of sizes, including an 80" x 48" version that fits perfectly in front of cart washers – #80524.

SURGICAL TRAY CORNER PROTECTORS

Prevent moisture retention and unwanted tears in sterile tray wraps with Surgical Tray Corner Protectors. Reticulated foam ensures maximum airflow and moisture management while providing a soft, supportive barrier to protect sterile wrappings. *(pictured to the left)*



"The single-use design reduced both our overall costs and our stress levels."

- Debbie Valdovinos, Sterile Processing Lead Tech
 Aurora St. Luke's in Milwaukee, WI



ISO 13485 CERTIFIED
 FDA REGISTERED

702 Prominence Drive
 New Kensington, PA 15068

Tel: 724-337-5500
 Fax: 724-337-0555

xodusmedical.com
 info@xodusmedical.com



LITBRCS003 R1 3/30/2020